

The Tube takes you to Pylorus Station


Calibration Tube for Sleeve Gastrectomy, Gastric Bypass and Mini Gastric Bypass



Temperature stable silicone tube

Technical Information:

- 1 soft conical radiopaque tip with suction hole
- 2 6 small holes avoid tissue suction damage
- 3 cm marking for precise placement
- 4 universal suction connector
- 5 clip to close tube for leak test
- 6 colour coded
- 7 vacuum control



TUBE32W	32 Fr
TUBE34R	34 Fr
TUBE36B	36 Fr
TUBE38G	38 Fr
TUBE40Y	40 Fr

