

The 1st choice in Calibration Tubes



Calibration



2019-01 EN

Technical Information:

- 1 non kinking, silicon tube
- 2 soft tip
- 3 cm marking
- 4 clip to close tube
- 5 hole in tip to avoid vacuum



REF	TUBE32W	32 Fr/Ch = 10,7 mm
REF	TUBE34R	34 Fr/Ch = 11,3 mm
REF	TUBE36B	36 Fr/Ch = 12,0 mm
REF	TUBE38G	38 Fr/Ch = 12,7 mm
REF	TUBE40Y	40 Fr/Ch = 13,3 mm



www.bariatric-solutions.com
info@bariatric-solutions.com

Phone: +41 - 52 741 - 65 00
Fax: +41 - 52 741 - 65 02