



The difference makes it happen.

MiniMizer Regular Gastric Band

MiniMIZER
REGULAR



 **BARIATRIC
SOLUTIONS**



First closing position 10.5cm (ø 31 mm)
 Second closing position 9.75cm (ø 26 mm)



- 1 two closing positions available:
less postoperative adjustments
- 2 inner reservoir integrated in gastric band
precise, predictable concentric filling
- 3 pre-shaped gastric band
for easy positioning
- 4 MMTIP Port
radiopaque, MRI Compatible

REF MMRE	MiniMizer Regular	Box 1 unit
REF MMTIP	MiniMizer Titanium Port Set	Box 1 unit

MiniMIZER

REGULAR



info@bariatric-solutions.com
 www.bariatric-solutions.com

Phone: +41 - 52 741 - 65 00
 Fax: +41 - 52 741 - 65 02