



The **Green**MAMBA retraction made easy.

- deployment in 4 turns
- validated disinfection & sterilisation process
- 8 cm (straight and with 45° angle)
and 12 cm wide
- non reflective shaft to prevent glares
- 5 mm diameter, 33 cm length



- REF **GMR8** 8 cm wide
- REF **GMR8B** 8 cm wide, 45° angle
- REF **GMR12** 12 cm wide
- REF **FCXL** FlushCYLINDER for cleaning

GMR-FL_p/Rev.02

CE

