

# A.M.I.<sup>®</sup>

## SOFT GASTRIC BAND PREMIUM SYSTEM



- comes in 2 sizes ( 35 and 32 mm diameter)
- 360° soft balloon
- safety closing mechanism
- low-pressure system
- unique catheter protection sleeve
- catheter separate from closing mechanism



### Port needle available

REF **AGB 382**

A.M.I. Port Needle

- Ø 20G / 0,9mm
- — 50<sup>mm</sup>

CE<sub>0297</sub>

REF **AGB 374\_LP**

Soft Gastric Band Premium Standard  
incl. Low Profile Port

REF **AGB 376\_LP**

Soft Gastric Band Premium Long  
incl. Low Profile Port

 **BARIATRIC  
SOLUTIONS**  
INTERNATIONAL



[www.bariatricsolutionsinternational.com](http://www.bariatricsolutionsinternational.com)  
[info@bariatricsolutionsinternational.com](mailto:info@bariatricsolutionsinternational.com)  
Phone: +41 52 741 15 59